

LOC440173-LOC440173 Gene

Reference: PTXBC027471

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized LOC440173		

For research use only.

GENE DESCRIPTION	LOC440173
GENE ID	440173
ORIGIN SPECIES	Human
ACCESSIONS	BC027471
VECTOR	pDONR223
NCBI SYMBOL	LOC440173
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaacccggtacacagatggaatgcagaaatcacccgtctctgcgtcgtcacgctggagctgtagaccggagcttattcgccatcttggctcctccctctatgagccaaccaggatggcgacttgactggattttacagttcattgacacgcatgagctctatcaaaattctg
ggagttttctgaaatgtttacatggaaggcacctggagccattga