

CNOT6-CCR4-NOT transcription complex, subunit 6 Gene

Reference: PTXBC027476

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CCR4; Ccr4a; CCR4-NOT transcription complex subunit 6; carbon catabolite repression 4 protein; carbon catabolite repressor protein 4 homolog; cytoplasmic deadenylase		

For research use only.

GENE DESCRIPTION	CCR4-NOT transcription complex, subunit 6
GENE ID	57472
ORIGIN SPECIES	Human
ACCESSIONS	BC027476
VECTOR	pDONR223
NCBI SYMBOL	CNOT6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgcatcatgacaattccagtgaaagtgagctggagctgggttgactaatgagactgaggagcagctttcctacgatctgcattatgtaacacaggtccagagagctttatggaagcgggagaggaggagcacttactcatgtgtattgttaaggaggatgtcatctttcatagatg
ctggaaactagagtgacattgttagatgctaaaggttgagctttacacaaaatgtctcatctgtattgttattgtctacaatatattga