

## PABPC3-poly(A) binding protein, cytoplasmic 3 Gene

Reference: PTXBC027617

|                |                                                                                                                                                                                                                                                                 |                     |                |
|----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                                                                                             | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                                                                                            | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | PABP3; PABPL3; tPABP; polyadenylate-binding protein 3; PABP-3; poly(A)-binding protein 3; poly(A)-binding protein-like 3; testis PABP; testis-specific poly(A)-binding protein (PABP); testis-specific poly(A)-binding protein 3; poly(A) binding protein cytop |                     |                |

For research use only.

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | poly(A) binding protein, cytoplasmic 3                                                           |
| <b>GENE ID</b>          | 5042                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC027617                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | PABPC3                                                                                           |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon)

atgaacccagcaccacccagctacccaacggcctcgtctacgtgggggacctccacccgacgtgactgagcgatgctctacgagaagttcagcccgccagggccctcctccatccggatctgcagggaactgacaccagcggtcctccaactacgcgtatgtgaactcca  
gcatacgaaggacgcgggacatgctctggacacatgaattttgatgtataaaggcgaagccagtcagcatcatgtggtctcagcgtgatccatctcgaaaaagtgagtgaggcaacataatctgtaaaatctggataagtcattataataaagcactgtatgatacagtttctgcttt  
ggtacacatccttctgttaacgtggtttgtgatgaaatggtccaagggttatggtattgtacatttgagacacacgaagcagctgaagagctattataaaaaatgaacggaatgctcctaaatggctcgcaaatgtttgtgacaatttaagtctctgaagaagcagaagcgaactgga  
gctaggcgcaaaagagttcccaatgtttacatcaagaattttgagagaacatggatgatgagcgcccttaaggatctcttggcgaagttcggcccgcccttaagtgtgaagtaatgaccgatgaaagtggaataatccaaggatttgatttgaagctttgaaagcatgaagatgcacaga  
aagctgtatgatgagatgaatgaaaggagctcaatggaaacaaattacgttggctcagctcagaaaaaagtggaacggcagacggaacttaagcgacattgaacagatgaagcaagataggatcaccagataccaggtgttaactttatgtgaaaaatcttgatgatggtattgatg  
atgaacgtctccgaaagcgttttctcatttggatcaatcactatgtgcaaaaggttatgaggaaggtggtcgcagcaaaaggtttgtttgtatgtttctcctccagaaagcactaaagcagttacagaaatgaacggtgaattgtggcccaaaagccattgtatgtagctttagctca  
gcgcaaaagaagcgcagcgttacctcactaacgagtatgacagaaatggcaagtgtagcagctgtgcccacacgcgagcacctcctcaggttactcatgacagctgtcccgacagctcagaacacatgctgcatactatcctcctagccaaattgctcgactaagaccaagtcctc  
ctggagctgctcaggggtgccagacctcatccattccaaaaaagccaggtgctatccgcccaggtgctcctagatgacatttagtactatgagacagctcttcacaggttccacgagtcagtcacgcagcgtgttgtaacacatcaacagacagtggttccacgtcctgcagctg  
ctgctgctgctgcagctacccctgctgtgcgcacgggttcacggtataaatatgctgcgggagttcgaatctcagcaacatcgtaatgcacagccacaagttacaatgaacagctgtgttcagtcataaaggtcaggaaacttgactgcctccaggttgccatctgccccctctcaaaa  
gcaaaagcaaatgttaggtgaacggctcttctcttattcaagccatgcaccctactctgctggtggaaaaactggtgagattgataattcagaactctttatgtctcagctccagagtcactccgttctaaggttgatgaagctgtagctgtactacaagccccaagctaa  
agaggctaccagaaagcagttacagctgctaccgggtgttccaactgtttaa