

MRPL23-mitochondrial ribosomal protein L23 Gene

Reference: PTXBC027710

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	L23MRP; RPL23; RPL23L; 39S ribosomal protein L23, mitochondrial; L23 mitochondrial-related protein; L23mt; MRP-L23; ribosomal protein related to L23 (mitochondrial); mitochondrial ribosomal protein L23		

For research use only.

GENE DESCRIPTION	mitochondrial ribosomal protein L23
GENE ID	6150
ORIGIN SPECIES	Human
ACCESSIONS	BC027710
VECTOR	pDONR223
NCBI SYMBOL	MRPL23
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgcgaatgtggtgtacccctgtaccgctgggtggccacaaacttcgggtgtccgaaccaactttcttcagctggtgcggccgggtgtggccagcccgaggacacgtgcagttccggatcccatggaatgacaagggtggacctcaggaattacctcagggcgc
tataacgtgccgtggctgctgtgctggacacgggtgcagcatggcttaacaagagaagatcacagaacgtgaggtatcaagaagccggactacaaggctgcctacgtgcagctggccatggacagaccttcacgttccagatctgtttccgagaaagacgagagccctgaag
gcagcgtgccgacgacctctacagcatgctcaggaggagaggcagagcagagcagcagcccgccggggcggtgtccctcagctggttcgggctgtga