



Reference: PTXBC027863

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SCR3; RNA-binding motif, single-stranded-interacting protein 2; suppressor of CDC2 with RNA-binding motif 3; RNA binding motif single stranded interacting protein 2		

For research use only.

GENE DESCRIPTION	RNA binding motif, single stranded interacting protein 2
GENE ID	5939
ORIGIN SPECIES	Human
ACCESSIONS	BC027863
VECTOR	pDONR223
NCBI SYMBOL	RBMS2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]