

HBZ-hemoglobin, zeta Gene

Reference: PTXBC027892

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HBZ-T1; HBZ1; hemoglobin subunit zeta; HBAZ; hemoglobin zeta chain; hemoglobin, zeta; zeta-globin		

For research use only.

GENE DESCRIPTION	hemoglobin, zeta
GENE ID	3050
ORIGIN SPECIES	Human
ACCESSIONS	BC027892
VECTOR	pDONR223
NCBI SYMBOL	HBZ
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctctgaccaagactgagaggaccatcattgtgtccatgtggccaagatctccacgcaggccgacaccatcggcaccgagactctggagaggtctctcctagccaccgcagaccaagactacttcccgaacttcgacctgcacccgggtccgcgcagttgcgcgcgcacgg
ctccaagtggtggccgccgtggcgacgcgggtgaagagcatcgacgacatcggcgccctgtccaagctgagcgcgctgcacgcctacatctgcgcgtggaccgggtcaactcaagctcctgtccactgcctgtgtcacctggccgcgcgttccccgcgacttcacg
gccgagggccacgcgcctgggacaaagtctatcgtgtctgtctgaccgagaagtaccgtga