

RHBDD1-rhomboid domain containing 1 Gene

Reference: PTXBC027900

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RRP4; rhomboid-related protein 4; rhomboid domain-containing protein 1; rhomboid-like protein 4; rhomboid domain containing 1		

For research use only.

GENE DESCRIPTION	rhomboid domain containing 1
GENE ID	84236
ORIGIN SPECIES	Human
ACCESSIONS	BC027900
VECTOR	pDONR223
NCBI SYMBOL	RHBDD1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcctgttttttttaaacgaaaggcgttttttccagtggtgttaccaggacggcaataactactttaatagttcaggcagctctggatcaggattattatccgcatggcagccagatcactatgaagaagcaccagggaactatgacacgtacacagcaggactgagtgaagaa
gaacagctcgagagagcattacaagccagcctctgggacggaggaataaccagaataagccaccaccctacgggtttcattctctaccagaagaatgaggagacagcggttcacagattcgatagccagtga