

PPP3R1-protein phosphatase 3 (formerly 2B), regulatory subunit B, alpha isoform Gene

Reference: PTXBC027913

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CALNB1; CNB; CNB1; calcineurin subunit B type 1; calcineurin B, type I (19kDa); protein phosphatase 2B regulatory subunit 1; protein phosphatase 2B regulatory subunit B alpha; protein phosphatase 3 (formerly 2B), regulatory subunit B (19kD), alpha isoform		

For research use only.

GENE DESCRIPTION	protein phosphatase 3 (formerly 2B), regulatory subunit B, alpha isoform
GENE ID	5534
ORIGIN SPECIES	Human
ACCESSIONS	BC027913
VECTOR	pDONR223
NCBI SYMBOL	PPP3R1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggaatgagcgaattatccttggaaatgtgctcacactttgatgcggatgaaatgaaaggctaggaaagagatttaagaagcttgattggacaattctggttcttgagtggtgaagagttcatgtctgcctgagttacaacagaatcctttagtacagcgagtaatatattcgac
acagatgggaatggagaagtagactttaagaattcattgagggcgtctcagttcagttcaaaaggagataaggagcagaattgaggtttgcttccgtatctgacatggataaagatggctatatttccaatggggaactctccaggtaattgaagatgatggtgggaacaatctgaa
agatacacagttacagcaattgtagacaaaccataataatgcagataaggatggagatggaagaataatccttgaagaattctgtgctgtttaggtggcctagatccacaaaagatgggtgtagatgtgtga