

DEFA3-defensin, alpha 3, neutrophil-specific Gene

Reference: PTXBC027917

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DEF3; HNP-3; HNP3; HP-3; HP3; neutrophil defensin 3; defensin 3, neutrophil-specific; defensin, alpha 3, neutrophil-specific; neutrophil peptide 3; defensin alpha 3		

For research use only.

GENE DESCRIPTION	defensin, alpha 3, neutrophil-specific
GENE ID	1668
ORIGIN SPECIES	Human
ACCESSIONS	BC027917
VECTOR	pDONR223
NCBI SYMBOL	DEFA3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaggaccctcgccatccttctgtgccattctctgtgtggccctgcaggccagctgagccactccaggcaagagctgatgaggttctgcagccccggagcagattgcagcggacatccagaagtgtgtttcccttgcattgggacgaaagcttgctccaaagcatccaggctcaaggaaaaacatggactgctattgcagaataccagcgtgcatfgcaggagaacgtcgctatggaacctgcactaccagggaagactctggcattctgctgctga