

CCL21-chemokine (C-C motif) ligand 21 Gene

Reference: PTXBC027918

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	6Ckine; CKb9; ECL; SCYA21; SLC; TCA4; C-C motif chemokine 21; Efficient Chemoattractant for Lymphocytes; beta chemokine exodus-2; chemokine (C-C motif) ligand 21; secondary lymphoid tissue chemokine; small inducible cytokine subfamily A (Cys-Cys), member		

For research use only.

GENE DESCRIPTION	chemokine (C-C motif) ligand 21
GENE ID	6366
ORIGIN SPECIES	Human
ACCESSIONS	BC027918
VECTOR	pDONR223
NCBI SYMBOL	CCL21
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctcagtcactggctctgagcctctatctgttctggttggccttggcctccagcagccaaggcagtgatggaggggtcagcagctgttgcctcaagtacagccaagggaagattcccgccaaggtgtccgcagctaccggaagcaggaaccaagcttaggctgtccatccagctatcctgttcttggccccgcaagcgctctcagcagagctatgtgcagacccaaggagctctgggtgcagcagctgatgcagcatctggacaagacaccatccccacagaaccagcccagggtgcaggaaggacagggggcctccaagactggcaagaaggaaagggtccaaaggctgcaagagactgagcgggtcacagaccctaaggccatag