

SNX20-sorting nexin 20 Gene

Reference: PTXBC027944

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SLIC1; sorting nexin-20; selectin ligand-interactor cytoplasmic 1; sorting nexin 20		
For research use only.			

GENE DESCRIPTION	sorting nexin 20
GENE ID	124460
ORIGIN SPECIES	Human
ACCESSIONS	BC027944
VECTOR	pDONR223
NCBI SYMBOL	SNX20
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcaagtcagagcaccctgggagccctggctgcatgggacccataaccagtcacggcgaaggaccagcaggaagaccagccactggccccgacctcccgacccagacactgacgggcaacttagacacacacagtgccctgagctccaactccagcatgaccacggg
gagcttcagcagtagctggcagaaccagaaatgccgtggaagcagtcacaaactgctctttgagatcgcttcagctcgcacgagagagaaaagtcttaagtgtgatggggaatcaaggcctggagagatgacttatccagggtcacgtgctgagacaggacagcaccagaacc
agaccgagatgtccagtcacaaagtacatgctctga