



Antibodies and Proteins

**ProteoGenix SAS**  
Espace Européen de l'Entreprise  
15 rue de La Haye  
67300 Schiltigheim, France  
(0)3 90 20 54 70 - Fax +33 (0)3 88  
contact@proteogenix.fr  
www.ProteoGenix.science

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88  
 contact@proteogenix.fr  
 www.ProteoGenix.science

## FCAR-Fc fragment of IgA, receptor for Gene

Reference: PTXBC027953

|                |                                                                                                                                           |                     |                |
|----------------|-------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                       | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                      | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | CTB-61M7.2; FcalphaRI; immunoglobulin alpha Fc receptor; Fc alpha receptor; Fc fragment of IgA, receptor for; Fc fragment of IgA receptor |                     |                |

**For research use only.**

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | Fc fragment of IgA, receptor for                                                                  |
| <b>GENE ID</b>          | 2204                                                                                              |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC027953                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | FCAR                                                                                              |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

**ORF (Only for reference, product provided without STOP codon )**

[illegible]

More info: Please contact us for any question: [contact@proteogenix.fr](mailto:contact@proteogenix.fr)