

CCL19-chemokine (C-C motif) ligand 19 Gene

Reference: PTXBC027968

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CKb11; ELC; MIP-3b; MIP3B; SCYA19; C-C motif chemokine 19; CC chemokine ligand 19; CK beta-11; EBI1-ligand chemokine; MIP-3-beta; beta chemokine exodus-3; chemokine (C-C motif) ligand 19; epstein-Barr virus-induced molecule 1 ligand chemokine; exodus-3; m		

For research use only.

GENE DESCRIPTION	chemokine (C-C motif) ligand 19
GENE ID	6363
ORIGIN SPECIES	Human
ACCESSIONS	BC027968
VECTOR	pDONR223
NCBI SYMBOL	CCL19
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccctgctactggccctcagctgctgttctctggactccccagcccccaactctgagtgaccacatgatgctgaagactgctgctgtgtgacccagaaacccctcgggtacatcgtgaggaactccactacctctcatcaaggatgctgcagggtgcctgctagtgtt
caccacactgagggccgccagctctgtgcacccccagaccagccctgggtagaacgcatcatccagagactgcagaggacctcagccaagatgaagcgccgcagcagttaa