



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

FKBP11-FK506 binding protein 11, 19 kDa Gene

Reference: PTXBC027973

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PPIase FKBP11; peptidyl-prolyl cis-trans isomerase FKBP11; 19 kDa FK506-binding protein; 19 kDa FKBP; FK506 binding protein 11, 19 kDa; FKBP-11; FKBP-19; rotamase; FK506 binding protein 11		

For research use only.

GENE DESCRIPTION	FK506 binding protein 11, 19 kDa
GENE ID	51303
ORIGIN SPECIES	Human
ACCESSIONS	BC027973
VECTOR	pDONR223
NCBI SYMBOL	FKBP11
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr