



Reference: PTXBC027985

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SAST; SAST170; microtubule-associated serine/threonine-protein kinase 1; syntrophin associated serine/threonine kinase; syntrophin-associated serine/threonine-protein kinase; microtubule associated serine/threonine kinase 1		

For research use only.

GENE DESCRIPTION	microtubule associated serine/threonine kinase 1
GENE ID	22983
ORIGIN SPECIES	Human
ACCESSIONS	BC027985
VECTOR	pDONR223
NCBI SYMBOL	MAST1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]