

C17orf45-chromosome 17 open reading frame 45 Gene

Reference: PTXBC027986

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C17orf45; NCRNA00188; TSAP19; C17orf76 antisense RNA 1		

For research use only.

GENE DESCRIPTION	chromosome 17 open reading frame 45
GENE ID	125144
ORIGIN SPECIES	Human
ACCESSIONS	BC027986
VECTOR	pDONR223
NCBI SYMBOL	C17orf45
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttccccggaagcctgtctagagggcgaggcagctgttgagatggcgtgctccccggtcctgcgcccgctggcgttctgcggcgccgtgacagcctggcaggcgagcgccggccgaatccggctgccgtggctgaggctcatggatcactctttt
gtgttagggccacatctgccagagcctggagctgcgaaggccgggacccgggtccccggccacagtgggggtgtgcaaacccgagagaactggattgcgtaccactgcagagtgctgaagacgggtagccacgaggttgcaaatctgtgaagaatcagcatcatgtttggcag
ctgagtattggagccaggagcctgcatgagggtttga