



Reference: PTXBC027990

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DRDL5813; PRO19630; tetratricopeptide repeat protein 14; TPR repeat protein 14; tetratricopeptide repeat domain 14		

For research use only.

GENE DESCRIPTION	tetratricopeptide repeat domain 14
GENE ID	151613
ORIGIN SPECIES	Human
ACCESSIONS	BC027990
VECTOR	pDONR223
NCBI SYMBOL	TTC14
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]