



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

C1orf51-chromosome 1 open reading frame 51 Gene

Reference: PTXBC027999

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

C1orf51; CHRONO; GM129; circadian-associated transcriptional repressor; chIP-derived repressor of network oscillator; computationally highlighted repressor of the network oscillator; circadian associated repressor of transcription

For research use only.

GENE DESCRIPTION	chromosome 1 open reading frame 51
GENE ID	148523
ORIGIN SPECIES	Human
ACCESSIONS	BC027999
VECTOR	pDONR223
NCBI SYMBOL	C1orf51
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr