

FLT3LG-fms-related tyrosine kinase 3 ligand Gene

Reference: PTXBC028001

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FLT3L; fms-related tyrosine kinase 3 ligand; flt3 ligand; fms related tyrosine kinase 3 ligand		

For research use only.

GENE DESCRIPTION	fms-related tyrosine kinase 3 ligand
GENE ID	2323
ORIGIN SPECIES	Human
ACCESSIONS	BC028001
VECTOR	pDONR223
NCBI SYMBOL	FLT3LG
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagcggctcaagactgtcgtgggtccaagatgcaaggcttgctggagcgcgtgaacacggagatacattgtcaccatgtgctttcagccccccagctgtttcgttcgttcagaccaacatctccgctcctgcaggagacctccgagcagctgggtggcgtgaagc
cctggatcactcgcagaactctccgggtgctggagctgcagtgctcagccgactcctcaacctgcccacccccatggagtcctccggccccctggagggccacagccccgacagccccgagccccctgtgctcctactgtgctgcccgtgggctcctgtgctgcccgtgctgct
gggtgctgactggcagaggacgcggcgaggacacccccctggggagcaggtgcccccgctcccgagtcctccagacctgctgctgtggagcaactga