



Reference: PTXBC028082

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD210; CD210a; CDW210A; HIL-10R; IL-10R1; IL10R; interleukin-10 receptor subunit alpha; IL-10 receptor subunit alpha; IL-10R subunit 1; IL-10R subunit alpha; IL-10RA; interleukin 10 receptor, alpha; interleukin-10 receptor alpha chain; interleukin-10 rece		

For research use only.

GENE DESCRIPTION	interleukin 10 receptor, alpha
GENE ID	3587
ORIGIN SPECIES	Human
ACCESSIONS	BC028082
VECTOR	pDONR223
NCBI SYMBOL	IL10RA
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]