

TRBV5-4-T cell receptor beta variable 5-4 Gene

Reference: PTXBC028083

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	T cell receptor beta variable 5-4
GENE ID	28611
ORIGIN SPECIES	Human
ACCESSIONS	BC028083
VECTOR	pDONR223
NCBI SYMBOL	TRBV5-4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccctgggctcctctgctgggcgtctgttctctgggagcaggctcagtgagactggagtcacccaagtcccacacactgatcaaacgagaggacagcaagtgacttgagatgctctctcagtcctggcacaacactgtctctgttaccacaggccctgggtcag
gggccccagttatctttagtattatagggaggaagagaatggcagaggaaactccctcctagattctcaggtctccagttccctaatatagctctgagctgaatgtgaacgcttggagctggagactcggccctgtatctctgtccagcagcatagcacgtgagcagttctcggcc
aggacacggctcaccgtgctagaggacctgaaaaacgtgttccacccgaggtcgtgtgtttgagccatcagaagcagagatctcccacccaaaaggccacactgtgtgtcctggccacaggcttctaccccgaccacgtggagctgagctggtggtgaatgggaaggagt
gcacagtgggtcagcagacacccgagccctcaaggagcagccgcccctcaatgactccagatactgctgagcagccgctgagggtctcggccacttctggcagaaccccgcaaccacttccgtgtcaagtcagttctacgggctctcggagaatgacaggtggacca
ggataggggccaaacctgtcaccagatctcagcgcggagccctgggttagagcagactgtggttcacctccagttaccagcaagggtctgtctgccaccctctatgagatcttctaggggaaggccacttctatgctgtgtgtcagtgccctgtgtgtgagccatg
gtcaagagaaggattccagaggctag