

BEST3-bestrophin 3 Gene

Reference: PTXBC028087

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	VMD2L3; bestrophin-3; vitelliform macular dystrophy 2-like 3; vitelliform macular dystrophy 2-like protein 3; bestrophin 3		

For research use only.

GENE DESCRIPTION	bestrophin 3
GENE ID	144453
ORIGIN SPECIES	Human
ACCESSIONS	BC028087
VECTOR	pDONR223
NCBI SYMBOL	BEST3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacaacagatgaaggaaattattcaaccacctaagctctcatctgaataattgggtccattcatctggttggaaatctgcaactaaagccggaatgaagtagaatcagagacagtgttgatctgcaatcattgatgactgaaatgaatcgataccgctcttggtgcagcctcttatt
cgggttatgactgggttgggattccgctggtttacacccaggtagcagagcagcttatcaacccttttggagaagatgatgatatttggaaactaactgtgcaattgacagaaatttgcaggtctctcttttagctgtggacgaaatgcacatgagcttaccgaatgaagaaggacattactgg
gacgattctgctgctgccaccatacacattggcagctgctgactactgcataccctcatttctgggtcaacagtcagatgggctgctggttccgacttctgacgaggagtggtgtgggattatgagaagcatggccatcggtccatccatgataagaagagtcagcggttctga
gtgcccacgaacacccctccagcccagaagaagctacaggaggcagacaagtgcagctcattgttcttaccctcagatgacctcagcccagccagggacactctggtatgcccctcaagaacccccccagggcctcaccacctggagaatcctgcttccagaaggaa
gccccacgctgcacttcagatgggagagctgtccaccatcagggagaccagccagacaagcactttatag