

SMEK1-SMEK homolog 1, suppressor of mek1 (Dictyostelium) Gene

Reference: PTXBC028088

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SMEK1; FLFL1; KIAA2010; MSTP033; PP4R3; PP4R3A; smk-1; smk1; serine/threonine-protein phosphatase 4 regulatory subunit 3A; SMEK homolog 1, suppressor of mek1; protein phosphatase 4 regulatory subunit 3A		

For research use only.

GENE DESCRIPTION	SMEK homolog 1, suppressor of mek1 (Dictyostelium)
GENE ID	55671
ORIGIN SPECIES	Human
ACCESSIONS	BC028088
VECTOR	pDONR223
NCBI SYMBOL	SMEK1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatgatgaagatgatgatgaggatgaagataaggaagatacggtaccattgtcaagaaagcaaaattgattcataataatggcaacggcctaggatcagctgttgaaaaaaactggtctccacccctccccatacaaaatccacaacaagcagtggtctctgtgaatgactga