

C6orf218-chromosome 6 open reading frame 218 Gene

Reference: PTXBC028118

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C6orf218; long intergenic non-protein coding RNA 518		
For research use only.			

GENE DESCRIPTION	chromosome 6 open reading frame 218
GENE ID	221718
ORIGIN SPECIES	Human
ACCESSIONS	BC028118
VECTOR	pDONR223
NCBI SYMBOL	C6orf218
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggattcagtgatgaacaggaagtcatcatgagttttggccacaccccaagtgacgacttagccagaatgggataactgggttccctactctctttatcatcctcaatgagagtaccaaattagagctagatggaaccttagtgaatactggctactgtccctccacca
gcctgccaccattcaagttgaagagacaaagacacatggaccttatgtaattactgggattacccaggaagtctgtggcaaaagtcagcttctccctccctgcttcccgccctgtctctgtgtacttttaccacactgggctgtttctgtatcacacttaa