

C3orf24-chromosome 3 open reading frame 24 Gene

Reference: PTXBC028122

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C3orf24; FANCD2 opposite strand protein; Fanconi anemia group D2 protein opposite strand transcript protein; FANCD2 opposite strand		

For research use only.

GENE DESCRIPTION	chromosome 3 open reading frame 24
GENE ID	115795
ORIGIN SPECIES	Human
ACCESSIONS	BC028122
VECTOR	pDONR223
NCBI SYMBOL	C3orf24
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcaggataccagctctgtgacccactggatgagatttccatggctgcgcacacgacacctacaccttcccaagcaccattcaaggctccccctgttccacacacacgctcgaccttgaagtgcagctgtgtttcaagggtcactctagctctagacag
ccccattcctggaatctggagtgtgagtcaccaagttaccctgccacacacagagttgcgcacgatgaacaacaaaggactgtcaggaagccccagccatccgctcagtgagtagattctgtctttggcaggggttatcacagctcagccacaaagtgaccgggactttcagagtttcag
acaagtcagccttttgcataatcattagcagggagcaccagtgcccatggactgaaggagcctcagattcagatgacagtcactatgtgcaaacagatgctgcgtctatcctcttgcgtgatgcaacttacaaaaagtcaccttgccttgagcactccaagtaa