



Reference: PTXBC028135

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HP2U; P2RU1; P2U; P2U1; P2UR; P2Y2; P2Y2R; P2Y purinoceptor 2; ATP receptor; P2U nucleotide receptor; P2U purinoceptor 1; P2U receptor 1; purinergic receptor P2Y, G-protein coupled, 2; purinoceptor P2Y2; purinergic receptor P2Y2		

For research use only.

GENE DESCRIPTION	purinergic receptor P2Y, G-protein coupled, 2
GENE ID	5029
ORIGIN SPECIES	Human
ACCESSIONS	BC028135
VECTOR	pDONR223
NCBI SYMBOL	P2RY2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]