

FLJ42289-hypothetical LOC388182 Gene

Reference: PTXBC028144

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical LOC388182
GENE ID	388182
ORIGIN SPECIES	Human
ACCESSIONS	BC028144
VECTOR	pDONR223
NCBI SYMBOL	FLJ42289
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacggccgcccgtttctggacagtgccattccagagccaggcaggagggcctcgggagggctggggctcccggtccgggtgaagcggagccgagcgcgccagaaccagaggaggaagcacacaggcccccagagagtggttgcaagagtcactcggggattag
gatccggccaaatcgcggaatttgaggctggaatcaagctcaacaggactgctgtgtgtaagcactcaagaggccctacaaccaacctgccaggagtcggctcctga