



Reference: PTXBC028158

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	OR2L14; olfactory receptor 2L13; olfactory receptor 2L14; olfactory receptor, family 2, subfamily L, member 14; olfactory receptor family 2 subfamily L member 13		

For research use only.

GENE DESCRIPTION	olfactory receptor, family 2, subfamily L, member 13
GENE ID	284521
ORIGIN SPECIES	Human
ACCESSIONS	BC028158
VECTOR	pDONR223
NCBI SYMBOL	OR2L13
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]