



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

C8orf38-chromosome 8 open reading frame 38 Gene

Reference: PTXBC028166

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UPF0551 protein C8orf38, mitochondrial; C8orf38; NADH dehydrogenase (ubiquinone) complex I, assembly factor 6; NADH:ubiquinone oxidoreductase complex assembly factor 6		

For research use only.

GENE DESCRIPTION	chromosome 8 open reading frame 38
GENE ID	137682
ORIGIN SPECIES	Human
ACCESSIONS	BC028166
VECTOR	pDONR223
NCBI SYMBOL	C8orf38
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgccgctccgcgcacggcctctctctggggcgccgttgaggcgttggcafcgcccgctgtgctgcgcggcgccgcctctgggtctgtacgcgcgcacggcgccgtgcccggggccggaggtgtctggggcgagcgtggctgcggccacgcggaccggcgccgtggggcactgac
cactactcgtcggagctgctgcggaaacgggattatgaaggttatctatctccctgctgctccctgcagaatccgaagctctgttttgacatgaggcccttaatgtggaactggtcaggctggtgattacttctgctgctgtctgtactgtatgccactgggactctgaacactaaactatgc
taa