

FOXR1-forkhead box R1 Gene

Reference: PTXBC028191

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DLNB13; forkhead box protein R1; forkhead box R1 variant 1; forkhead box R1 variant 2; forkhead box protein N5; forkhead box R1		

For research use only.

GENE DESCRIPTION	forkhead box R1
GENE ID	283150
ORIGIN SPECIES	Human
ACCESSIONS	BC028191
VECTOR	pDONR223
NCBI SYMBOL	FOXR1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggggaacgagctctttctggcctcaccacatcacctccccttagcggagcagaaactgccaggtataaactccgaattgttaagccacaaaattacccctagagaaaaaaccaaccctgataaggatgggtccagattatgagcccaacctctggatgtggtaaatcccaattg
tgatccccctggaagctggaggtctcaggacgttaggaagaggaggacctgacaagcacactcccctctctcagccaccccaagagggaagatgccagctgctcagaggccgaggggtggaatcactgtcccagtcctccagcaagcgggtctccccctggaagcgggttg
ccctttccccagcacctgggagctcacagaagaggaggaggtctgaggaccaggaagacagctctctatggctctccatccccacaaaaggccccctccagagtcggaggtctggcaagccagcagcaggcggggaggtctgtgtcccgccccctcaattactcca
cctaattgccctggcattaaagaacagttccccctgtggcctcaactgcaacagatctacagtttctcctgaaagcacttcccccttttccggacggccccggaaggctggaagaatactgtccgtcacaatctctgttttcgagacagcttggagaaagtgcctgtcagcatcaggcgagg
gccagcacacggcctcgtatcttgcctctggaagttgaccgaggaggacacccgcttgcggaggagggccgcttgccttcactcggtagaagatccaacagtgcatgagccagcagatgtgatgcccttctcttgcatttaa