



Reference: PTXBC028207

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PER; RIGUI; hPER; period circadian protein homolog 1; Period, drosophila, homolog of; circadian clock protein PERIOD 1; circadian pacemaker protein RIGUI; hPER1; period homolog 1; period circadian clock 1		

For research use only.

GENE DESCRIPTION	period homolog 1 (Drosophila)
GENE ID	5187
ORIGIN SPECIES	Human
ACCESSIONS	BC028207
VECTOR	pDONR223
NCBI SYMBOL	PER1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr