



Reference: PTXBC028212

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	F730038I15Rik; FOAP-2; P101-PI3K; p101; phosphoinositide 3-kinase regulatory subunit 5; PI3-kinase p101 subunit; phosphatidylinositol 4,5-bisphosphate 3-kinase regulatory subunit; protein FOAP-2; ptdIns-3-kinase p101; phosphoinositide-3-kinase regulatory		

GENE DESCRIPTION	phosphoinositide-3-kinase, regulatory subunit 5
GENE ID	23533
ORIGIN SPECIES	Human
ACCESSIONS	BC028212
VECTOR	pDONR223
NCBI SYMBOL	PIK3R5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]