

ZNF471-zinc finger protein 471 Gene

Reference: PTXBC028227

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ERP1; Z1971; zinc finger protein 471; EZFIT-related protein 1		
For research use only.			

GENE DESCRIPTION	zinc finger protein 471
GENE ID	57573
ORIGIN SPECIES	Human
ACCESSIONS	BC028227
VECTOR	pDONR223
NCBI SYMBOL	ZNF471
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaggagaaccagtcgcaaaagcctgaaattccaaaaccagaatgcctctctcctccaaaggatcacactcctgccggcaagggaacaaaactggatagagaatgagttgacaaattgacagaagtaggctcagaagggtgggaaataacaaactcctctgagctaaaggag
catgttctaaccgaatgcaaggagcctaagaacctgaaaaagggttagagggaattgctaactagataaccagtttagagaagaacataaatgacctgatggagctgaaaaacacagcatga