



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

C18orf56-chromosome 18 open reading frame 56 Gene

Reference: PTXBC028301

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C18orf56; TYMS opposite strand protein; TYMS opposite strand		
For research use only.			

GENE DESCRIPTION	chromosome 18 open reading frame 56
GENE ID	494514
ORIGIN SPECIES	Human
ACCESSIONS	BC028301
VECTOR	pDONR223
NCBI SYMBOL	C18orf56
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagcccccctcgggggcccagcgcatactggggcgactgggagcccgcccgaggagccctggggagcggcttaacctccggctgctgctgtgtgtgtgtggcccgccagtaacgtccctaattggaccctcgtttcggcgtctgaagcagaggaggacatgctccctccctggcaggaatgaggttaggactaaacggggctccgagcgcgccggcagctcccgagcgtctccccgcgcgcctcctctctccgcacacccgtccgcaggggcccgggcgtaacctctcagcgtgagcgcgcctgatgccgaataaccacagggtgcgcgtgccg