

GFER-growth factor, augmenter of liver regeneration Gene

Reference: PTXBC028348

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ALR; ERV1; HERV1; HPO; HPO1; HPO2; HSS; FAD-linked sulfhydryl oxidase ALR; ERV1 homolog; erv1-like growth factor; hepatic regenerative stimulation substance; hepatopoietin protein; growth factor, augmenter of liver regeneration		

For research use only.

GENE DESCRIPTION	growth factor, augmenter of liver regeneration
GENE ID	2671
ORIGIN SPECIES	Human
ACCESSIONS	BC028348
VECTOR	pDONR223
NCBI SYMBOL	GFER
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcggacgcagcagaagcgggacaccaagttagggaggactgccgcccggatcgcgagggaactggcgccacacagctgggtgtcctccacacctggccgctactacccgacctgcccacccagaacagcagcaagacatggccagttcatacatttttaagttta
ccccctgtgaggagtgctgaagacctaaagaaaggctgtgcagggaaccacccagacacccgacccggcgatgcttcacacagtggtgtgcccactgcacaaatgaagtgaaccgcaagctggcgcaagcctgacttcgactgctcaaaagtggatgagcgtggcgcgacgctg
gaaggatggctcctgtgactag