

SLC37A3-solute carrier family 37 (glycerol-3-phosphate transporter), member 3 Gene
Reference: PTXBC028380

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	sugar phosphate exchanger 3; solute carrier family 37 (glycerol-3-phosphate transporter), member 3; solute carrier family 37 member 3		

For research use only.

GENE DESCRIPTION	solute carrier family 37 (glycerol-3-phosphate transporter), member 3
GENE ID	84255