



Reference: PTXBC028385

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	rotamase PPIL2; CYC4; CYP60; Cyp-60; UBOX7; hCyp-60; peptidyl-prolyl cis-trans isomerase-like 2; PPIase; U-box domain containing 7; cyclophilin, 60kDa; cyclophilin-60; cyclophilin-like protein CyP-60; peptidylprolyl cis-trans isomerase; peptidylprolyl iso		

GENE DESCRIPTION	peptidylprolyl isomerase (cyclophilin)-like 2
GENE ID	23759
ORIGIN SPECIES	Human
ACCESSIONS	BC028385
VECTOR	pDONR223
NCBI SYMBOL	PPIL2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1