

LIN28-lin-28 homolog (C. elegans) Gene

Reference: PTXBC028566

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	LIN28; CSDD1; LIN-28; ZCCHC1; lin-28A; protein lin-28 homolog A; RNA-binding protein LIN-28; zinc finger CCHC domain-containing protein 1; zinc finger, CCHC domain containing 1; lin-28 homolog A		

For research use only.

GENE DESCRIPTION	lin-28 homolog (C. elegans)
GENE ID	79727
ORIGIN SPECIES	Human
ACCESSIONS	BC028566
VECTOR	pDONR223
NCBI SYMBOL	LIN28
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggtcctgttccaaccagcagtttgcaggtggctgcgccaaggcggcagaagagggcggcgaggagggcgaggagacggcgccggcgggcgaggagcctcagctgtgcacgggtgcgggcatctgtaagtggttcaacgtgcgcatgggggttcggttcctgtccatga
ccgcccgcgcgggggtcgcgctcgaacccagtgatgtcttgcacacagagtaagctgcacatggaagggttcggagcttgaaggaggggtgagggcagtgagggttcacctttaagaagtcagccaagggtctggaatccatccgtgcacgggacctgggtgagatctgtattg
ggagtgagagggcgccaaaaggaaagagcatgcagaagcgcagatcaaaaggagacaggtgtctacaactgtggaggtctagatcatcatgccaaggaatgcaagctgccacccagcccaagaagtgccacttctccagagcatcagccatggttagcctcatgtcgcctgaag
gcccagcaggccctagtcacaggggaaagccaacctcttcgagaggaagaagaataccacagccctaccctgtcccgaggcacagaattga