

C1orf21-chromosome 1 open reading frame 21 Gene

Reference: PTXBC028567

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C1orf21; PIG13; cell proliferation-inducing gene 13 protein; proliferation-inducing protein 13; chromosome 1 open reading frame 21		

For research use only.

GENE DESCRIPTION	chromosome 1 open reading frame 21
GENE ID	81563
ORIGIN SPECIES	Human
ACCESSIONS	BC028567
VECTOR	pDONR223
NCBI SYMBOL	C1orf21
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggctgtgctcgcgaagcatgttccactgttcaaatgaagaggaagcccagaagggaactaccagaacggagatgtttggcgatgagtataatgatcaaacccagtggaagaggtcaatacatgaaaaatgggcagaagaagagcagaaatagcagccaggaac
caagaaacttgaaaaagtgccagctcaatgtaagacttaaaactaataagaggttcgggattagttcatcaaccagagcaaacatgcacatctctgaaagccaacaagaattctcagaatgctgatgaaaaatgaaaagggtcgggattactgttcggaagaagagatat
cacatag