



Reference: PTXBC028607

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TDE1L; TDE2; TMS-2; TMS2; serine incorporator 1; tumor differentially expressed protein 2		

For research use only.

GENE DESCRIPTION	serine incorporator 1
GENE ID	57515
ORIGIN SPECIES	Human
ACCESSIONS	BC028607
VECTOR	pDONR223
NCBI SYMBOL	SERINC1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]