

ZCCHC24-zinc finger, CCHC domain containing 24 Gene

Reference: PTXBC028617

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C10orf56; Z3CXXC8; zinc finger CCHC domain-containing protein 24; zinc finger, 3CxxC-type 8; zinc finger, CCHC domain containing 24; zinc finger CCHC-type containing 24		

For research use only.

GENE DESCRIPTION	zinc finger, CCHC domain containing 24
GENE ID	219654
ORIGIN SPECIES	Human
ACCESSIONS	BC028617
VECTOR	pDONR223
NCBI SYMBOL	ZCCHC24
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagcctgtgtcggccatcgacacgagcgccgctcgtgtaccagcccgccagctgctcaactgggtctacgtgtcgtcgcaggacacgcaccaggtctagcgccttcgatgccttcggccgtgccgaccgcccggcgccgcacccccggagctggccttggcaagggccg
ccccgagcagctgggtcgcctcgcctcgcactccagctatctcaacagcttctccagctgcagcgaggagggcgtgagcaacagtggtgtacaagggcgccctcacctatggctccctcaacaacatcgccgatggcctcagctccctcaccgagcacttctcagacctgacctcacctc
cgaggtcgcgaagccagcagcgcccccacccaactacgtgtccacctgtgttcaacaaggacactacatcaaggactgccccagggcagccccaaaggcgagggcctgactccataccaggcaaaaagcgctgtctcggcgagtaaatgttcccaagtgcagagaa
aatggatgagcgggaactcctgggccaacatggggcaggagtgcatcaagtgcacatcaactgtgtatccacacaagcagagaccctggagaagcccagcgctggagtggtccgaccagagcagcagcccgagcagcactctgcgagaagtgcagggtcctgggtacta
ctgccgtcgcgtgcagtga