



Reference: PTXBC028682

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HSD8; p56; coiled-coil domain-containing protein 91; GGA-binding partner; p56 accessory protein; coiled-coil domain containing 91		

For research use only.

GENE DESCRIPTION	coiled-coil domain containing 91
GENE ID	55297
ORIGIN SPECIES	Human
ACCESSIONS	BC028682
VECTOR	pDONR223
NCBI SYMBOL	CCDC91
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]