



Reference: PTXBC028717

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ARNO3; GRP1; PSCD3; cytohesin-3; ARF nucleotide-binding site opener 3; PH, SEC7 and coiled-coil domain-containing protein 3; general receptor of phosphoinositides 1; pleckstrin homology, Sec7 and coiled-coil domains 3; cytohesin 3		

For research use only.

GENE DESCRIPTION	cytohesin 3
GENE ID	9265
ORIGIN SPECIES	Human
ACCESSIONS	BC028717
VECTOR	pDONR223
NCBI SYMBOL	CYTH3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]