



Reference: PTXBC028974

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ZNF424; zinc finger protein 57; zinc finger protein 424		
For research use only.			

GENE DESCRIPTION	zinc finger protein 57
GENE ID	126295
ORIGIN SPECIES	Human
ACCESSIONS	BC028974
VECTOR	pDONR223
NCBI SYMBOL	ZNF57
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]