



Reference: PTXBC029041

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HOF; ODA-8S; PRIMS; ZNF288; zinc finger and BTB domain-containing protein 20; BTB/POZ zinc finger protein DPZF; dendritic cell-derived BTB/POZ zinc finger; zinc finger protein 288; zinc finger and BTB domain containing 20		

For research use only.

GENE DESCRIPTION	zinc finger and BTB domain containing 20
GENE ID	26137
ORIGIN SPECIES	Human
ACCESSIONS	BC029041
VECTOR	pDONR223
NCBI SYMBOL	ZBTB20
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]