

STC1-stanniocalcin 1 Gene

Reference: PTXBC029044

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	STC; stanniocalcin-1; stanniocalcin 1		
For research use only.			

GENE DESCRIPTION	stanniocalcin 1
GENE ID	6781
ORIGIN SPECIES	Human
ACCESSIONS	BC029044
VECTOR	pDONR223
NCBI SYMBOL	STC1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtccaaaactcagcagtgcttctggtgctgtgatcagtgcttctgcaaccatgagcgagcagaatgactctgtgagccccaggaaatcccgagtgccgctcaaaactcagctgaagtgggtgctgctcaacagtgctctacaggtcggctggggctttgcatgcctgga
aaactccacctgtgacacagatgggatgtatgacatctgtaaatctctgttacagcgtctgctaaattgacactcagggaaaagcattcgtcaagagagcttaaaatgcatcgccaacggggtcacctccaaggtctctcgcctcgcattcggaggtgctccactttccaaagatgattgtg
aggtgcaggaagagtgctacagcaagctgaatgtgtgagcatcgccaagcggaaacctgaagccatcactgaggtcgtccagctgcccgaatcacttctccaacagatactataacagactgtccgaagcctgctggaatgtgatgaagacacagtcagcacaatcagagacagcctga
tggagaaaattggcctaactggccagccttccacatcctgcagacagaccactgtgccaacacacccacgagctgacttcaacaggagacgcaccaatgagccgcagaagctgaaagtctctcaggaacctccgaggtgaggaggactctccctccacatcaaacgcac
atcccatgagatgcataa