



Reference: PTXBC029125

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ALXR; FMLP-R-II; FMLPX; FPR2A; FPRH1; FPRH2; FPRL1; HM63; LXA4R; N-formyl peptide receptor 2; FMLP-R-I; FMLP-related receptor I; LXA4 receptor; RFP; formyl peptide receptor-like 1; lipoxin A4 receptor (formyl peptide receptor related); formyl peptide rece		

GENE DESCRIPTION	formyl peptide receptor 2
GENE ID	2358
ORIGIN SPECIES	Human
ACCESSIONS	BC029125
VECTOR	pDONR223
NCBI SYMBOL	FPR2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]