

C6orf35-chromosome 6 open reading frame 35 Gene

Reference: PTXBC029130

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UPF0463 transmembrane protein C6orf35; C6orf35; BM033; transmembrane protein 242		

For research use only.

GENE DESCRIPTION	chromosome 6 open reading frame 35
GENE ID	729515
ORIGIN SPECIES	Human
ACCESSIONS	BC029130
VECTOR	pDONR223
NCBI SYMBOL	C6orf35
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagacagcggcgctgcaactggcgagccgctctgggctggaggtccggggtccacgaatgacggcctttctgttaaagggtgaatttctgtaccgttgctgcagcgggaatgctagctggattattacaacattatcattggctaaaaagaaagccctgaatggctc
aataagggaagtatggccacggctgcattaccgaaagcgggtcttccctgcttgcgagctctgggctgggctccctgtatgcatgggtgtggggttgattagcttcgagctctggaagcttaggagttcacagatgaacgactttcgaagtaaatgcaatcaatattccaac
aattccaagaactccgaatcggctgtgagtgagggaacattgaaatccaatga