

PRR15-proline rich 15 Gene

Reference: PTXBC029131

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	proline-rich protein 15; proline rich 15		
For research use only.			

GENE DESCRIPTION	proline rich 15
GENE ID	222171
ORIGIN SPECIES	Human
ACCESSIONS	BC029131
VECTOR	pDONR223
NCBI SYMBOL	PRR15
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccgacagcggcgatgctggcagctccggccctgggtggaatcgtcaccaacagcagaagaaaaagcaaggaagccgagtggtgggtgcccctcccagcccgctccggggagcccacgcacctgcgcgccagcccgactggaccagcagctccggg
agaaccagcaccacaatctctcggggcgccggcgagcccccaaacagacaagtatacggggacaaatccggcagcagccgcgcaattgaagatctcgcgctccggcgctttaaggagaagagaaagtgcgcgccacgctgctccggaggcggcaggtcccg
aggaggcaggctttctggtgaccccccacgaggacaagcagtag