

CCDC42-coiled-coil domain containing 42 Gene

Reference: PTXBC029224

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CCDC42A; coiled-coil domain-containing protein 42A; coiled-coil domain containing 42		

For research use only.

GENE DESCRIPTION	coiled-coil domain containing 42
GENE ID	146849
ORIGIN SPECIES	Human
ACCESSIONS	BC029224
VECTOR	pDONR223
NCBI SYMBOL	CCDC42
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagtcctggcgcacatggaagaggaagacctggccgagtagtccggctgcagtagtggggagcggctgctgcagatgctccagaactcccaatgttgagggggcgtcggagtcgccatccatctggctactggagaagaaaaaggagacagaaatcatgcatcaaatatggtgca
gaagaagaagatgttcagcgcagaatggaacccctgaacctgcctgggaggaactggcgcttaaggaaagcccaactgaaggctcacatccagaagcttgagcagttcatccaggagaacgaccagaacggatccgcgcctgaagaagccaacaggagcgagaactcaag
tgccagcacatgcaggagctgaccaagcgcaagcaggagatggtggcgctgcggctggagcaccagcggctgagcgccaagctgaaggactactacatcttcaacaagtacctagagaaggtggtggagaactccgaggaatctcgtggcgccacatccagaacaccgcagcca
agaagaccctcctgctggcaccattaagatggccacgtgaacctcttcagatcgtgagcaagcacctgaaggaggtgactgaggtggcactggagacacccacaagcagctggacatgatccagcaatttatccaagaccggtcggacatctggcagaggtgaaaaagaagga
acaacagcggagtcggattaa